



iProven Releases Accurate FMH 139 5pack H Cg Pregnancy Home Test

October 07, 2019

October 07, 2019 - PRESSADVANTAGE -

iProven, a business based in Dordrecht, The Netherlands, has announced that they have released the FMH-139 5pack kit to test for pregnancy that they have ensured to have optimal sensitivity. With this new midstream test kit for pregnancy, those who are trying to get pregnant can determine with confidence whether they are finally pregnant.

This kit is designed to be used on cycle day one, which is the first day of the next menstruation cycle of the woman who is trying to conceive. When used on that particular day, the test kit is expected to provide results that are over 99 percent accurate. The FMH-139 5pack kit is provided with 5 test strips, which means that the tester can double check, triple check the results, or even more.

The company had developed this highly accurate midstream test kit for pregnancy so that women who are TTC may not suffer the intense disappointment from a false positive result. Imagine what a big let down it can be to feel so elated when the home test kit provides positive results, only to find out it should have been negative when they go to their doctor.

It would take 14 days or two weeks from ovulation to cycle day one, or from the date of the baby dance that resulted in fertilization up to the time when women TTC can apply the home pregnancy test. If the test is not accurate and it happens to provide a positive result, the disappointment could be quite severe.

Suzanne, a spokesperson for iProven, says, "We are so happy to have finally released the FMH-139 Midstream Pregnancy Test kit. It gives you the reassurance you need. From the day of your expected period, the pregnancy test will give you results over 99% accurate. The pregnancy test has optimal sensitivity to detect the pregnancy hormone hCG in your urine. With a sensitivity of 25 mIU/ML, you can rely on the accuracy of the iProven FMH-139 pregnancy test."

Meanwhile, the FMH-139 pregnancy test kit is provided with a pink cap. When doing the test, this pink cap has to be transferred to the other end. In this way, there is more clearance between the urine stream and the hand holding the test kit at one end. This allows the user to avoid that messy situation where she gets urine in her hand. Also, when performing the test, the user should wash her hands to ensure hygiene. Then, she should direct the urine stream onto the test strip for 10 seconds. After that, she simply waits for three to five minutes.

When taking a look at the results, the user should ensure that there is a control line, which is located on the right side of the window. If this is not present, the results would be invalid. If there is a control line but nothing more, the results would be negative. But if there are two lines, made up of the control line and the test line, it is time to celebrate because the results are positive. Moreover, the tester is over 99 percent sure that the results are correct.

It should be noted that it takes two weeks or 14 days from the date of conception for the body to produce sufficient amounts of the hCG pregnancy hormone for this to be detected in the urine. That is why it is recommended to use the test kit on the cycle day one. This is because after the egg has been fertilized, it has to travel through the fallopian tube into the uterus where it will attach itself to the uterine wall. After the fertilized egg has attached itself to the uterine wall, the body starts to produce a large amount of hCG. That is the time when the pregnancy hormone can be detected in the urine.

iProven is a company that develops home health products that are designed for various members of the family. People who would like to know more about the FMH-139 hCG pregnancy home test kit or other products offered by iProven can check out the company's website or their Amazon page. Or they can contact them on the telephone or through email.

###

For more information about iProven, contact the company

here:iProvenSuzanne1-503-974-0913marketing@iproven.comiProvenBeaverton, OR

iProven

It is our passion to develop high-quality products that can easily be operated by every home user and provide the best-tailored education together with the products. We love to help you get the most out of the product and be more health aware.

Website: <http://www.iproven.com>

Email: marketing@iproven.com

Phone: 1-503-974-0913



Powered by PressAdvantage.com