



iProven Digital Forehead Non Contact Thermometer Back In Stock

October 01, 2019

October 01, 2019 - PRESSADVANTAGE -

iProven, a business based in Dordrecht, The Netherlands, has announced that the NCT-336 BL non contact digital thermometer for babies, kids, and adults is now in stock after being temporarily unavailable for some time. Those who have been looking for something like the NCT-336 BL digital thermometer that can be used to measure body temperature without actually touching the person may have been disappointed when it became unavailable. Fortunately, it is now back in stock and people, especially parents, can have a way of checking if their babies or young kids have a fever without actually touching them and possibly irritating them.

“This thermometer has a sturdy design and gives you reliable results, over and over again. A thermometer you can count on!” says Suzanne, a spokesperson for iProven.

She continues, “And what is really convenient is that you don’t have to touch the forehead of your baby while measuring. Keep a distance of 2-6 inches between the thermometer and your baby’s forehead. Your little one won’t have time to be fussy, and you know how fussy babies and little tots are, particularly if they are not feeling well.”

The NCT-336 BL medical thermometer is easy to use. The user simply points it at the forehead of a baby or child at a distance of two to six inches to determine the body temperature within just two seconds. Furthermore, iProven has ensured that it is sturdy and reliable so that people can always rely on the readings from this thermometer. And they can take reading after reading and be sure it will always function properly, that is, provide accurate temperature readings. It even has a 100 days money back guarantee, which means that customers who are not satisfied with it can simply return it within 100 days after purchase, no questions asked.

What is more is that the NCT-336 BL non-contact thermometer can be used to measure the temperature of just about anything. Parents who are worried that the feeding bottle is too hot can check its temperature by simply putting the thermometer to Object Mode and then pointing it at the bottle from a distance of two to six inches. When worried about the baby's bath being too hot or too cold, parents can simply point it at the bath from a distance of two to six inches, click on it and the temperature of the bath is provided. It can also be used to measure room temperature, for parents who are worried that their babies might find it uncomfortable.

Also, the NCT-336 BL is convenient to use, particularly for babies and young kids. Parents know how irritable little ones can get if they are sick or not feeling well. But because there is no need to touch the young ones with the thermometer and because it can provide an accurate reading within two seconds, parents struggling to take the temperature of fussy babies and little tots are a thing of the past. And if the little one is asleep, mommies and daddies can simply mute the thermometer so that it will not emit a sound. And because it is non-contact, babies and little kids who are sleeping won't be disturbed in their sleep. Parents can simply take note of the green display and go back to sleep.

Another important feature is that the color of the display will change into yellow if the temperature is elevated and into red if there is a fever. There will also be an alarm, except when the mute function is on, so that parents are immediately alerted if there is a fever so that they can take the appropriate action.

iProven develops home health products just like the NCT-336 BL that can be used by people to help in caring for their own health and that of their family members. Those who require more information about the digital forehead thermometer and other similar products from iProven can view their website or their Amazon page. Those interested may also want to contact them on the telephone or via email.

###

For more information about iProven, contact the company here:iProvenSuzanne1-503-974-0913marketing@iproven.comiProvenBeaverton, OR

iProven

It is our passion to develop high-quality products that can easily be operated by every home user and provide the best-tailored education together with the products. We love to help you get the most out of the product and be more health aware.

Website: <http://www.iproven.com>

Email: marketing@iproven.com

Phone: 1-503-974-0913



Powered by PressAdvantage.com